


1. Part Number: 675 8051H1100S

2. Photograph:



3. Product Information:

Description:	Jumbo Pipette 1100mm - sterile Single use, disposable liquid sampler
Use:	
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none">• Conforms to FDA 21 CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
	A blue dye is added to one of the tips. The blue dye is certified: <ul style="list-style-type: none">• FDA acceptable for food contact• Free of materials of animal origin
Method of Construction:	Ultrasonic welding
Nominal Length	1090mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	Handle End (Blue) = 1.5mm Product End (Natural) = 3mm
Nominal Weight of Sampler:	123g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilization:	Gamma Irradiation (25kGy minimum dose)
Number of Samplers per Box:	37
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested