


1. Part Number: 675 8051H 450

2. Photograph:



3. Product Information:

Description: Jumbo Pipette 450mm Single use,
disposable liquid sampler High Density
Use:
Material of Construction: Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

A blue dye is added to the Tip. The blue dye certified as being FDA acceptable and does not contain any materials of animal origin.

Method of Construction: Ultrasonic welding
Nominal Length: 420mm
Diameter of Outer Tube: 21mm
Diameter of hole in tip: Handle End (Blue) = 1.5mm
Product End (Natural) = 3mm
Nominal Weight of Sampler: 50g
Moulding Environment: Class 100,000 Cleanroom
Assembly & Primary Packing Environment: Class 10,000 Cleanroom
Individually Bagged? Yes (heat sealed PE/PET bag)
Method of Sterilization: N/A
Number of Samplers per Box: 50
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This product has not been ATEX tested