



1. Part Number: 675 8051H450S

2. Photograph:



3. Product Information:

Description: Jumbo Pipette 450mm - Sterile Single use,

Use: disposable liquid sampler High Density

Material of Construction: Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

A blue dye is added to the Tip. The blue dye certified as being FDA acceptable and does not contain any materials of animal origin.

Method of Construction: Ultrasonic welding

Nominal Length 420mm

Diameter of Outer Tube: 21mm

Diameter of hole in tip: Handle End (Blue) = 1.5mm
Product End (Natural) = 3mm

Nominal Weight of Sampler: 50g

Moulding Environment: Class 100,000 Cleanroom

Assembly & Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilization: Gamma Irradiation (25 kGy minimum dose)

Number of Samplers per Box: 50

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This product has not been ATEX tested