

1. Part Number: 675 8063H1000S

2. Photograph:



3. Product Information:

| | |
|--|--|
| Description: | Powder Lance - 1000mm - Sterile Single |
| Use: | use, disposable powder sampler HDPE (High |
| Material of Construction: | Density Polyethylene) - virgin <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free |
| Method of Construction: | Ultrasonic welding |
| Nominal Length: | 953mm |
| Blade Length: | 851mm |
| Diameter of Blade: | 22mm |
| Nominal Weight of Sampler: | 120g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Welding & Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Quantity per Box: | 20 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This device has not been ATEX tested |