

1. Part Number: 675 8063H300S

2. Photograph:



3. Product Information:

Description:	Powder Lance - 400mm – Sterile Single use,
Use:	disposable powder sampler HDPE (High
Material of Construction:	Density Polyethylene) - virgin
	<ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	401mm
Blade Length:	284mm
Diameter of Blade:	22mm
Nominal Weight of Sampler:	64g
Moulding Environment:	Class 100,000 Cleanroom
Welding & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested