

1. Part Number: 675 8063H500S

2. Photograph:



3. Product Information:

Description:	Powder Lance - 500mm - Sterile Single
Use:	use, disposable powder sampler HDPE
Material of Construction:	(High Density Polyethylene) - virgin <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	601mm
Blade Length:	500mm
Diameter of Blade:	22mm
Nominal Weight of Sampler:	80g
Moulding Environment:	Class 100,000 Cleanroom
Welding & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 – 45kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested