

# Specification Sheet

Specification Number:

Revision:

10

Date Issued:

11-Jul-24

**1. Part Number:** 675 8060H1001

**2. Photograph:**



### 3. Product Information:

Description:	Powder Thief - 1000mm – 1ml Tip Single
Use:	use, disposable powder sampler HDPE
Material of Construction:	(High Density Polyethylene) - virgin <ul style="list-style-type: none"><li>• Conforms to FDA CFR 177.1520</li><li>• Conforms to EU Regulations 10/2011</li><li>• Conforms to EC Regulations 1935/2004</li><li>• BSE/TSE free</li></ul>
Method of Construction:	Ultrasonic welding
Nominal Length:	987mm
Diameter of Outer Tube:	21mm
Nominal Tip Size:	1ml
Width of Handle:	103mm
Nominal Weight of Sampler:	165g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	20
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested