

1. Part Number: 675 8060H501

2. Photograph:



3. Product Information:

Description:	PowderThief - 500mm – 1ml Tip Single
Use:	use, disposable powder sampler HDPE
Material of Construction:	(High Density Polyethylene) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	511mm
Diameter of Outer Tube:	21mm
Nominal Tip Size:	1ml
Width of Handle:	103mm
Length 'L':	5mm
Nominal Weight of Sampler:	93g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	20
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested