

1. Part Number:

675 8060H505

2. Photograph:



3. Product Information:

Description:	PowderThief - 500mm – 5ml Tip – Sterile Single
Use:	use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	511mm
Diameter of Outer Tube:	21mm
Nominal Tip Size:	5ml
Width of Handle:	103mm
Length 'L':	31mm
Nominal Weight of Sampler:	93g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	NA
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested