



1. Part Number: 675 8073H1000

2. Photograph:



3. Product Information:

Description: Slot Sampler - 1000mm Single use,
 Use: disposable powder sampler High Density
 Material of Construction: Polyethylene (HDPE) - virgin
 • Conforms to FDA 21 CFR 177.1520
 • Conforms to EU Regulations 10/2011
 • Conforms to EC Regulations 1935/2004
 • BSE/TSE free
 Method of Construction: Ultrasonic welding
 Overall Length: 968mm
 Diameter of Outer Tube: 26mm
 Diameter of Inner Blade: 22mm
 Nominal Weight of Sampler: 234g
 Slot Positions:



Number of Sampling Slots: 3
 Sampling Slot Width: 8mm
 Assembly & Primary Packing Environment: Class 100,000 Cleanroom
 Individually Bagged? Yes (heat sealed PE/PET bag)
 Method of Sterilization: N/A
 Number of Samplers per Box: 20
 Recommended Storage Conditions: Dry and ambient temperature
 Shelf Life: 5 years from the date of manufacture
 ATEX: Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application