



1. Part Number: 675 8073H300

2. Photograph:



3. Product Information:

Description:	Slot Sampler - 300mm
Use:	Single use, disposable powder sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Overall Length:	420mm
Diameter of Outer Tube:	26mm
Diameter of Inner Blade:	22mm
Nominal Weight of Sampler:	100g
Slot Positions:	



Number of Sampling Slots:	3
Sampling Slot Width:	11mm
Assembly & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture